

SIRIUS

60W THULIUM FIBER LASER




STONE & BPH

- EXTREME FREQUENCY EMISSION
- LOW RETROPULSION
- EFFECTIVE LITHOTRIPSY
- GREAT CUTTING & ABLATION CAPABILITY






SIRIUS



STONE SOLUTION

BPH SOLUTION



Reduced Size of Dust



Low Retropulsion



Extreme Frequency Emission



Compact Design



Reduced noise levels



Standart Electrical Installation



BPH



Hemostasis




Soft Tissues




FIBERS

Single Use & Reusable




FOOT PEDALS

INTUITIVE SCREEN



SIRIUS

SOLUTION



Technical Informations


Laser type	Thulium Fiber Laser
Wavelength	1940 nm \pm 20nm
Power of Target	60 Watts
Energy/Pulse	up to 6J
Application mode	Pulsed and CW
Time function	Pulsed Mode and CW
Pulse Duration	75 μ s < laser pulsion duration < 1000 μ s
Repetition rate	Up to 2500 Hz
Peak power	Up to 500W


Ordering Informations


Reference	Description
MS01	SIRIUS 60w Thulium Medical Laser
FIBERS	Single Use & Reusable

This brochure is not intended for the U.S. market.
Certain Intended Uses/Configuration/Models/Accessories are not cleared for U.S.

For more information visit www.rocamed.com

 Rocamed SAM
9 avenue Albert II
98000 Monaco

 (+377) 97 98 42 43

 info@rocamed.com

 www.rocamed.com