


Articulation lever  
Switch between self-lock  
mode and non-self-lock  
mode



**OTU<sup>®</sup> Medical**

Silicon Valley Based Medical Technology Company

## CONTACT


OTU Medical Inc.

2847A Whipple Rd., Union City, CA 94587, USA

[www.otumed.com](http://www.otumed.com)

[sales@otumed.com](mailto:sales@otumed.com)

+1 (408) 461-9801




# WiScope<sup>®</sup> Single-Use Digital Flexible Ureteroscope

8.6Fr

7.5Fr

6.3Fr

Bring smart views to One Time  
Use (OTU) surgical devices



**OTU<sup>®</sup> Medical**

- No cross-infection
- No sterilization
- No maintenance




Bring smart views to One Time Use (OTU) surgical devices

Model	OTU-110S OTU-110R	OTU-130S OTU-130R	OTU-140S OTU-140R
Field of View	100°		
Direction of View	Forward		
Depth of field	2 – 50 mm		
Insertion Tube Outer Diameter	8.6Fr	7.5 Fr	6.3 Fr
Working Channel Diameter	3.6 Fr	3.6 Fr	3.3 Fr
Angulation Range	275° Up / 275° Down		
Working Length of Shaft	670mm		

**OTU**® Medical

## WiScope® Single-Use Digital Flexible Ureteroscope



**6.3<sub>Fr</sub> | 7.5<sub>Fr</sub> | 8.6<sub>Fr</sub>**  
Insertion Tube Outer Diameter

**400×400**  
pixels

Up to **275°** bidirectional  
non-planar bending