







BoneScalpel Tip Selection Guide



Tip	Part #	Approach	Purpose	Procedures
 10mm Blade, Blunt .5 x 5 mm	MXB-10 Working Length* = 34 mm	Open	Osteotomy En-Bloc Resection Bone Cutting	Laminoplasty Laminectomy Facetectomy PSO TLIF/PLIF Corpectomy Laminotomy Foraminotomy ACDF
	MXB-10LS Working Length* = 145 mm	MIS		
	MXB-10LC Working Length* = 145 mm	MIS		
 20mm Blade, Blunt 1.0 x 6.4 mm	MXB-20 Working Length* = 40 mm	Open		
	MXB-20LC Working Length* = 145 mm	MIS		
	MXB-MIS-20LCRS Working Length* = 100 mm	Mini-Open/MIS		
 25mm Blade, Blunt 1.0 x 6.4 mm	MXB-25 Working Length* = 45 mm	Open		
 20mm Blade, Serrated 1.0 x 6.4 mm	MXB-B1 Working Length* = 40 mm	Open		
 Micro Hook Shaver 1.3 x 1.8 mm	MXB-S1 Working Length* = 57 mm	Open		
 Diamond Shaver 4.4 mm	MXB-S3 Working Length* = 34 mm	Open	360° Bone Shaving & Bone Shaping Decompression Gross Resection	Corpectomy VCR Decortication Laminotomy

*Approximate distance from typical index finger placement to end of probe using BoneScalpel Handpeice(BCM-HP)

If you would like further information or would like to evaluate BoneScalpel, please contact us at +1.631.694.9555



1938 NEW HWY, FARMINGDALE, N.Y. 11735 | +1.631.694.9555 +1.631.694.3285 FAX
 MISONIX, INC. | NASDAQ SYMBOL. MSON | MISONIX.COM

