

# TraxerFlow™ II

## Irrigation Line



EndoFlow™ II  
Mono & Double Chamber






» » » UROLOGY




## TraxerFlow™ Dual Port

- ➔ **PRESSURE CHAMBER**
  - Anti-reflux valve prevent back flow of fluid.
- ➔ **MANUAL PUMP**
  - Provide extra flow for clearing the surgical field.
- ➔ **LUER LOCK CONNECTION**
  - Grips the devices and prevent leakage.
- ➔ **REMOVABLE TUBING EXTENSION**
  - Allow surgeon to use the pump itself or with the assistant.
- ➔ **TRANSPARENT BLUE TUBE**
  - Allow to see the saline water.



### GRAVITY LINE

Référence	Product description
ONGT2000ST12	TraxerFlow DUAL PORT Irrigation Line for Ureteroscopy Laser



Compatible with all principal brands of endoscopes

### WITH ENDOFLOW II

Référence	Product description	
SOUR9107ST12	TraxerFlow Irrigation Line for Ureteroscopy Laser	1 spike
SOUR9207ST12	TraxerFlow Irrigation Line for Ureteroscopy Laser	2 spikes



MKG-FY-TFLOW-EN-00  
 Designed by FOCALYS. Illustrations in this document are non-contractual.  
 ROCAMED reserves the right to change product specifications without notice.

For more information visit [www.rocamed.com](http://www.rocamed.com)

Rocamed SAM  
 9 avenue Albert II  
 98000 Monaco

(+377) 97 98 42 43  
 (+377) 92 05 61 50 (fax)

info@rocamed.com

[www.rocamed.com](http://www.rocamed.com)