

Rocamed Sam
9, avenue Albert II
Le Copori
98000 Monaco

Tel : +377 97 98 41 99
Email : Info@rocamed.com

www.rocamed.com

✓ Urology

Technical Sheeth

Ref. (s)	
Intended use	
Description	
Main Materials	
Pack factor / Packaging	
Sterilization	
Manufacturer	
Regulations	
Shelf Life	



RocaTub Ureteral Catheter

ROTAxxxST - ROTGxxxST - ROTHxxxST- ROTPxxxST

RocaTub catheters are intended for use during ureteral catheterisation, drainage, opacification of the upper urinary tract (retrograde ureteropyelography) and flushing procedures.

The RocaTub catheters are sterile and for single use only. Please refer to labeling for the specifications of each product.

ROCAMED's RocaTub ureteral catheter series offers a wide range of catheters with various distal tip shapes.

Ureteral Catheter in PU, 80cm Length, radiopaque, graduated each 50mm, Supplied with a stylet and luer..

Distal Tip : Flexible Straight/Open // No Eye				
ROTP5000ST				
Distal Tip : Straight/Open // 2 Eyes				
4 CH	5 CH (4,8)	6CH	7CH	8CH
ROTA4200ST	ROTA5200ST	ROTA6200ST	ROTA7200ST	ROTA8200ST
Distal Tip : Straight/Open // No Eye				
4 CH	5 CH (4,8)	6CH	7CH	8CH
ROTG4000ST	ROTG5000ST	ROTG6000ST	ROTG7000ST	ROTG8000ST
Distal Tip : Curved/Open // No Eye				
4 CH	5 CH (4,8)	6CH	7CH	8CH
ROTH4000ST	ROTH5000ST	ROTH6000ST	ROTH7000ST	ROTH8000ST

Main Materials	PU	
Pack factor / Packaging	ROTxxxxST	Sterile peel pouch (45 x 140mm) 10 units per box (890 x 78 x 38mm) - 6 boxes per brown box (905 x 175 x 140mm)
Sterilization	Sterilized by Ethylene Oxide	
Manufacturer	ROCAMED France - Z.I de Signes - Allée de Stockholm - 83870 Signes - France	
Regulations	Certified ISO13485:2016 In compliance with Directive 93/42 / EEC MDD Class: IIa	
Shelf Life	3 years	

